TAUREON

NOVACOL®

FIBRILLAR

NL	Absorbeerbaar Hemostatisch Collageen	GEBRUIKSAANWIJZING	4
DE	Absorbierbarer Kollagenhämostaseverband	GEBRAUCHSANWEISUNG	7
UK	Absorbable Collagen Hemostat	DIRECTIONS FOR USE	10
FR	Collagène hémostatique résorbable	NOTICE D'UTILISATION	13
IT	Emostatico al collagene assorbibile	ISTRUZIONI PER L'USO	16
ESP	Colágeno hemostático absorbible	INSTRUCCIONES DE USO	19
GR	Απορροφήσιμο αιμοστατικό κολλαγόνο	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	22
TU	Emici Hemostatik Kollajen	KULLANMA TALİMATI	25

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DIRECTIONS FOR USE

DESCRIPTION OF THE PRODUCT

Novacol Fibrillar Absorbable Haemostatic Collagen is a purified bovine collagen product.

The material is sterile, pyrogen-free and absorbable. The haemostat action largely depends on the product's fibrous structure. This allows a rapid initial packing of the wound, which mechanically stops the bleeding.

Stabilisation of the initial packing is provided by the formation of a fibrin network. Once collagen comes into contact with blood, the platelets stick together and release coagulation factors that, together with plasma factors, lead to fibrin formation and ultimately to the formation of a clot.

INDICATIONS

Novacol Fibrillar can be used during surgical procedures, except in urologic or ophthalmologic surgery, as a haemostasis aid if the control of bleeding by a ligature or other conventional method is not having the desired effect or is unusable.

CONTRA-INDICATIONS

Novacol Fibrillar should not be used when closing skin incisions, because it can interfere with the skin's healing process. This interference can be ascribed to a normal mechanical interposition of dry collagen and not to some intrinsic interference with wound healing.

Novacol Fibrillar should not be applied to bone surfaces to which prosthesis materials are being attached with methylmethacrylate. The collagen in another absorbable collagen product has been reported to reduce the adhesive force of methylmethacrylate by occluding the pores in porous bone.

SPECIAL WARNINGS

Novacol Fibrillar is inactivated by autoclaving. It must not be resterilised. Just like with any other foreign substance, it can lead to the development of infections in contaminated wounds.

Novacol Fibrillar should not be used in situations involving pulsing arterial bleeding.

Novacol Fibrillar should not be used at sites where blood or other fluids are combined or in situations where the site of the bleeding is submerged. Novacol Fibrillar does not act as

a tampon or plug on the bleeding site, nor does Novacol Fibrillar close off the area where blood has collected behind a tampon.

Only the quantity of Novacol Fibrillar required for achieving haemostasis should be used. The long-term consequences of leaving Novacol Fibrillar in situ are unknown. Once opened, unused Novacol Fibrillar should be destroyed, because it is not possible to once again sterilise it.

Just as with other haemostatic products, it is not advisable to leave Novacol Fibrillar in an infected or contaminated space, nor should it be used in people known to be allergic to products of bovine origin. If Novacol Fibrillar is placed in cavities or closed spaces, then large quantities of Novacol Fibrillar should not be inserted, as it can absorb fluids, expand and exert pressure on the surrounding structures.

The safety of this product has not yet been determined in children and pregnant women. As a result, Novacol Fibrillar should only be used in those cases if the benefit of its application clearly exceeds the risks.

Novacol Fibrillar should not be used to treat systemic bleeding disorders. Novacol Fibrillar can disintegrate when used in water or a saline solution, thus producing small particles that can flow through the filters of equipment used for blood recovery. The same disintegration can also occur during an operation when applying the fibrous material to a bleeding site. If a blood recovery apparatus is being used in such a situation, then these small collagen particles can pass through the filter in the machine and end up in the patient's bloodstream when the recovered blood is transfused back. This can result in intravascular thrombus formation. Therefore, Novacol Fibrillar should not be used in sites where the use of a blood recovery apparatus is being considered.

ADVERSE REACTIONS

Novacol Fibrillar is a collagen product.

Reported adverse effects of other collagen haemostatic products include haematomas, infections, wound dehiscence, inflammation and oedema. Other reported adverse effects that might be connected to the use of collagen haemostatic products include adhesion formation, allergic reactions, foreign body reactions and, in a single case, subgaleal seroma. The use of microfibre collagen in dental cavities has been reported to improve the prevention of alveolalgia. It cannot be excluded that all of the above-mentioned reactions could arise when using Novacol Fibrillar.

APPLICATION

Novacol Fibrillar can be moulded (by hand or using an instrument) into the shape required by the surgical site and is directly applied to the bleeding surface while exerting pressure. Both the amount of pressure and the length of time it should be exerted depend on the severity and the nature of the bleeding to be stopped.

The time required to stop bleeding is determined by the type of surgical procedure and the severity of the bleeding before treatment. Usually, the bleeding is stopped with Novacol Fibrillar in 2 to 6 minutes.

If necessary, Novacol Fibrillar can be left in situ. The surgeon is responsible for deciding whether to remove surplus Novacol Fibrillar before closing the wound.

PACKAGING AND CONTENTS

Novacol Fibrillar is supplied in a non-moulded form in plastic containers in the following format:

0.5 g	IN BOXES OF 5	ARTICLE NO. 2010
1.0 g	IN BOXES OF 5	ARTICLE NO. 2011

The sterility of this product is guaranteed, unless the individual packaging is damaged or opened.

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