# TAUREON

## NOVACOL®

## PAD

NL	Absorbeerbaar Hemostatisch Collageen	GEBRUIKSAANWIJZING	4
DE	Absorbierbarer Kollagenhämostaseverband	GEBRAUCHSANWEISUNG	7
UK	Absorbable Collagen Hemostat	DIRECTIONS FOR USE	10
FR	Collagène hémostatique résorbable	NOTICE D'UTILISATION	13
IT	Emostatico al collagene assorbibile	ISTRUZIONI PER L'USO	16
ESP	Colágeno hemostático absorbible	INSTRUCCIONES DE USO	19
GR	Απορροφήσιμο αιμοστατικό κολλαγόνο	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	22
TU	Emici Hemostatik Kollajen	KULLANMA TALİMATI	25

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#### **DIRECTIONS FOR USE**

#### **DESCRIPTION OF THE PRODUCT**

Novacol Pad Absorbable Haemostatic Collagen is a purified bovine collagen product.

The material, prepared in the form of a non-woven bandage, is sterile, pyrogen-free and absorbable. The haemostatic action largely depends on the product's fibrous structure. This allows a rapid initial packing of the wound, which mechanically stops the bleeding. The formation of a fibrin network provides the stabilisation of the initial packing. Once collagen comes into contact with blood, the platelets stick together and release coagulation factors that, together with plasma factors, lead to fibrin formation and, ultimately, to the formation of a clot.

#### **INDICATIONS**

Novacol Pad can be used during surgical procedures (except in urologic or ophthalmologic surgery) as a haemostasis aid if the control of bleeding by a ligature or other conventional method is not having the desired effect or is unusable. In neurological procedures, such as laminectomy and craniotomy, the Novacol Pad can be applied to extradural sites only if the dura is intact or where direct contact with neural tissue can be avoided. The safety of this product with subdural structures has not yet been ascertained.

#### **CONTRA-INDICATIONS**

Novacol Pad should not be used when closing skin incisions, because it can interfere with the skin's healing process. This interference can be ascribed to a normal mechanical interposition of dry collagen and not to some intrinsic interference with wound healing.

Novacol Pad should not be applied to bone surfaces to which prosthesis materials are being attached with methylmethacrylate. The collagen in another absorbable collagen product has been reported to reduce the adhesive force of methylmethacrylate by occluding the pores in porous bone.

#### **SPECIAL WARNINGS**

Novacol Pad is inactivated by autoclaving. It must not be resterilised. Just as with any other foreign substance, it can lead to the development of infections in contaminated wounds.

Novacol Pad should not be used in situations involving pulsing arterial bleeding.

Novacol Pad should not be used at sites where blood or other fluids are combined or in situations where the site of the bleeding is submerged. Novacol does not act as a tampon or plug on the bleeding site, nor does Novacol close off the area where blood has collected

behind a tampon.

Only the quantity of Novacol Pad required for achieving haemostasis should be used. The long-term consequences of leaving Novacol in situ are unknown. Once opened, unused Novacol Pads should be destroyed, because it is not possible to once again sterilise them.

Just as with other haemostatic products, it is not advisable to leave a Novacol Pad in an infected or contaminated space, nor should it be used in people known to be allergic to products of bovine origin. If a Novacol Pad is placed in cavities or closed spaces, then large quantities of Novacol Pad should not be inserted, as it can absorb fluids, expand and exert pressure on the surrounding structures.

The safety of this product has not yet been determined in children and pregnant women. As a result, Novacol Pad should only be used in those cases if the benefit of its application clearly exceeds the risks.

Novacol Pad should not be used to treat systemic bleeding disorders.

## **ADVERSE REACTIONS**

Novacol Pad is a collagen product. Although various postoperative complications have been noted in patients treated with Novacol, the investigator could ascribe none of these complications to Novacol.

Reported adverse effects of other collagen haemostatic products include haematomas, infections, wound dehiscence, inflammation and oedema. Other reported adverse effects that might be connected to the use of collagen haemostatic products include adhesion formation, allergic reactions, foreign body reactions and subgaleal seroma (in a single case). The use of microfibre collagen in dental cavities has been reported to improve the prevention of alveolalgia. It cannot be excluded that all of the above-mentioned reactions could occur when using Novacol Pad.

## **APPLICATION**

Novacol Pad is directly placed on the bleeding surface, while exerting pressure.

Novacol Pad can be cut to measure. Both the amount of pressure and the length of time it should be exerted depend on the severity and the nature of the bleeding to be stopped.

The time required to stop bleeding is determined by the type of surgical procedure and the severity of the bleeding before treatment. Usually, the bleeding is stopped with Novacol Pad in 2 to 6 minutes.

Novacol Pad retains its structure in the presence of blood and does not disintegrate in a moist environment. It can be simply removed from the site once the bleeding has stopped. Novacol Pad works best when applied dry.

If necessary, Novacol Pad can be left in situ. The surgeon is responsible for deciding whether to remove surplus Novacol Pad before closing the wound. Animal studies have shown that the absorption and tissue reactions to Novacol Pad are identical to those reactions reported with another absorbable haemostatic collagen product.

### CLINICAL STUDIES

The safety, effectiveness and characteristic ease of use of Novacol Pad have been extensively evaluated in a range of surgical procedures.

The average time to achieve haemostasis with Novacol was 4.4 minutes. The Passive Haemagglutination Assay (PHA) and Enzyme-Linked Immuno-absorbent Assay (ELISA) methods were used to assess the immunological potential of Novacol to stimulate antibody production in animals. These analyses showed little or no evaluation of antibody-titres in animals treated with Novacol, while a slight increase in antibodies was noted in animals treated with a control collagen haemostatic agent. The ease of use of the Novacol Pad was considered superior to that of the control collagen haemostatic agent.

### PACKAGING AND CONTENTS

Novacol Pad is supplied in a non-woven form in plastic film packages in the following formats:

2.5 cm x 5 cm	IN BOXES OF 10	ARTICLE NO. 3001
5 cm x 8 cm	IN BOXES OF 10	ARTICLE NO. 3002
8 cm x 10 cm	IN BOXES OF 10	ARTICLE NO. 3003

The sterility of this product is guaranteed, unless the individual packaging is damaged or opened.

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